




Product Overview




NumaLink™

Solves multi-vendor compatibility problems. Seamlessly transfer patient studies between all nuclear medicine and PET systems or to PACS. Acquire on one vendors camera or workstation and process and display on another.




NumaList™

NumaList connects to the HIS/RIS and adds or corrects DICOM modality worklist information in images. Used with NumaLink, files can be translated into DICOM with the corrected worklist prior to sending to PACS.




NumaList Plus™

NumaList Plus enhances NumaList by adding full DICOM Worklist capabilities to the Philips/ADAC Pegasys®. This eliminates time-consuming re-entry of patient information.




NumaStore™

Expandable, secure, reliable, image management and archive solution, optimized for NM, PET/CT and preclinical imaging departments. This Vendor Neutral Archive automates workflow throughout the department with support for DICOM, non-DICOM and preclinical image formats. Available as a complete system or software only.



NumaRead™

Migrate offline tape/MOD/CD image archives to NumaStore or PACS. Supports multiple media types and storage devices. Automatically moves the entire media contents to another system or provides quick and easy access to images on tape or disk. Works with NumaLink to translate non-DICOM images to DICOM. Rental options are available.



Custom Configurations

Numa has the expertise to optimize the most complex imaging department's workflow. Numa's sales and engineering team will examine your department's requirements and respond with a time and cost saving proposal. Multiple vendors, imaging systems, locations or a mix of older and new technologies are Numa's specialty. Contact Numa today!



For more information please contact Numa Sales at 603-883-1909 or sales@numa-inc.com

