

IHE Integration Statement			
Vendor	Product Name	Version	Date
Numa, Inc.	NumaStatus		1/19/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Consistent Time	Time Client	None
Radiation Exposure Monitoring for Nuclear Medicine	Radiopharmaceutical Activity Supplier	Radiopharmaceutical Activity Supplier Worklist Query
Radiation Exposure Monitoring for Nuclear Medicine	Radiopharmaceutical Activity Supplier	None
Radiation Exposure Monitoring for Nuclear Medicine	Radiopharmaceutical Activity Supplier	Patient Based Worklist Query

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	http://www.numa-inc.com/resources/numastatus_dicom_conformance_statement.pdf

Link to general information on IHE : www.ihe.net