

NumaBridge Host System Requirements:

The NumaBridge software requires the following minimum host system specifications.

OPERATING SYSTEM	CPU		RAM	HDD Space
	# Cores	Speed		
Windows 7, 10 Pro (32 or 64bit) Windows Server 2008 R2, 2012 R2, 2016	2	2.0ghz	4GB	1GB

- Two (2) Network Interface Card(s) (Two NICs are required to isolate Windows 2000/XP systems from the network)
- Keyboard, display (1024x 768 or better) and Mouse; compatible KVM switch or remote control application such as RDP. NumaBridge does not require the user to access the user interface.
- Internet/VPN access for Numa support.